SWISS MADE 

# SERIES LT

Brief instructions / Kurzanleitung / Courte notice d'utilisation / Instrucciones breves

Precisa Gravimetrics AG


**Precisa Gravimetrics AG**  
**Moosmattstrasse 32**  
**CH-8953 Dietikon**  
**Tel.: +41-44-744 28 28**  
**Fax.: +41-44-744 28 38**  
**Email: [service@precisa.ch](mailto:service@precisa.ch)**

<http://www.precisa.com>

### Customer service

Information on and addresses for local customer service centres can be found on our website.

### Copyright

These operating instructions are protected by copyright. All rights reserved. No part of these operating instructions may be reproduced, edited, copied or distributed in any form either through photocopying, microfiche, reprinting or any other method, and in particular via any electronic means, without written permission.

© Precisa Gravimetrics AG, 8953 Dietikon, Switzerland, 2014.

	<b>Table of contents</b>		
<b>1.</b>	<b>Introduction .....</b>	<b>6</b>	
<b>2.</b>	<b>Safety precautions .....</b>	<b>6</b>	
<b>3.</b>	<b>Technical data .....</b>	<b>7</b>	
<b>4.</b>	<b>Starting up .....</b>	<b>8</b>	
4.1	Choice of a suitable location ..	8	
4.2	Unpacking the balance .....	8	
4.3	Assembling the balance.....	9	
4.4	Under-floor weighing .....	10	
4.5	Levelling the balance.....	11	
4.6	Calibrating the balance .....	11	
<b>5.</b>	<b>Key functions.....</b>	<b>12</b>	
<b>6.</b>	<b>Overview of menus .....</b>	<b>13</b>	
6.1	Menu control .....	13	
6.2	Configuration menu.....	14	
6.3	Application menu.....	15	
<b>7.</b>	<b>Key menus .....</b>	<b>16</b>	
<b>7.1</b>	<b>Tare menu .....</b>	<b>16</b>	
7.1.1	Taring the balance.....	16	
7.1.2	Other functions .....	16	
<b>7.2</b>	<b>Mode menu .....</b>	<b>17</b>	
7.2.1	Switch to an activated application or back to weighing mode .....	17	
7.2.2	Other functions .....	17	
<b>7.3</b>	<b>Print menu .....</b>	<b>18</b>	
7.3.1	Print individual value / report.....	18	
7.3.2	Other functions .....	18	

**DANGER**

Warning of a possible danger which can lead to death or serious physical injuries.

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**WARNING**

Warning of a possibly dangerous situation which can lead to less severe injuries or damage.

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**NOTE**

Tips and important rules on the correct operation of the balance.

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


## 1.

### Introduction

There are various models in the 321 Series. These brief instructions describe the balance with windshield. If you are using a less complex model, the descriptions of elements that are not included can be skipped. These brief instructions are intended as an initial overview and are an extract from the original operating instructions. The brief instructions must be available at all times at the workstation. More detailed descriptions can be found in the original operating instructions, which can be downloaded from our website. You should read the original operating instructions in full to ensure that you make the best possible use of your balance.

## 2.

### Safety precautions

- When the balance is used in environments with increased safety requirements, the relevant precautions must be taken.
- The balance must only be operated with a fault-free power supply unit.
- Before connecting the power supply unit, ensure that the operating voltage specified on the power supply unit matches the mains supply.
- The balance must only be operated with the power supply unit supplied.
- If safe operation of the balance is not possible, it must be disconnected from the power supply and secured against unintentional use.

## 3.

**Technical data**

Casing dimensions W x D x H for: - Readability 0.1 mg - Readability 1 mg - Readability > 1 mg	228 x 381 x 345 mm 228 x 381 x 124 mm 228 x 381 x 90 mm
Net weight: - Readability 0.1 mg - Readability 1 mg - Readability > 1 mg	6.6 kg 5.5 kg 5.6 kg
Weighing range for each model (see details on nameplate)	120 g to 10 200 g
Allowable ambient conditions (for calibratable balances, see details on the nameplate)	5 °C to 40 °C
Allowable relative humidity	25 % - 85 %, non-condensing
Mains connection (19)	90 V - 264 V, 47 Hz - 63 Hz
Power consumption	18 W, 12 V DC, 1.5 A
Data transfer connection (17 and 18)	RS 232/V24 interface, USB type B (device)

Tab. 2  
Technical data

## 4. Starting up

### 4.1 Choice of a suitable location

- Solid, firm and low-vibration, horizontal base
- No direct sunlight
- No draft
- Ambient temperature: 5 °C - 40 °C
- Relative humidity (non-condensing): 25 % - 85 %

### 4.2 Unpacking the balance

#### **NOTE**

Retain the original packaging in order to avoid damaging the balance during shipping or transport and to enable the balance to be stored in ideal conditions when out of use for long periods.

- Carefully unpack the balance.


## 4.3

## Assembling the balance

**CAUTION**

Before connecting the device to the mains, be sure that the operating voltage stated on the device and on the power adaptor agrees with the local mains voltage. Do not touch any key during the following start-up sequence until after about 3 seconds an acoustic signal sounds.


- Place the dust protection cover on the balance.
- Install the draft shield (1).
- Place the balance pan holder in position and add the weighing pan (2).
- Insert the supplied power supply unit into the socket (9) on the back of the balance.

## 4.4 Under-floor weighing



For weighing objects which, because of their size or shape, cannot be put in or on the balance.

- Remove the screw (10) on the underside of the balance.
- Place the balance over an opening.
- Hang the hook for under-bench weighing, to which the object can be secured, into the aperture (12) of the metal casting (11).


## 4.5 Levelling the balance

- Adjust the two adjustable feet (6) so that the air bubble in the level control (5) below the buttons on the control panel (3) lies precisely in the centre of the sight glass marking.

Wrong

Right

## 4.6 Calibrating the balance

- Press the  key to switch on the balance.
- Press the  key until the 20-character info line (14) in the display (4) shows «CALIBRATION».

The balance begins the calibration process.

- Follow the instructions in the info line (14) and the eight-digit measurement display (13).

### NOTE

The balance must be recalibrated each time its location is changed.

## 5. Key functions

Key	Designation	Function	Position index
	ON/OFF	Switches the balance on/off	15
	MENU / esc	Calling up the configuration menu and the application menu	16
	O/T / ins O/T / enter	Activates the tare functions and calibration function and other functions	17, 20
	CHANGE / clr	Switches between the basic program and the chosen application	19
	PRINT	Activates printing functions and the Back button in the menu	18
	FUNCTION KEYS left / right up / down	Function keys. Start the functions in the info-line	21, 22, 23, 24

Tab. 5  
Key functions

## 6.

## Overview of menus

**NOTE**

With user guidance activated, the keys that currently have a function or can be used backlit.

## 6.1

**Menu control**

Navigating through menu levels / Making an input:

Key	Designation	Function
 	left / right	Changes from menu main path into sub-paths and back
 	up /down	Up/Down movements within main and sub-paths Changes selected parameters
	O/T / enter	Selects parameters Stores the changed parameters
	MENU / esc	Interrupts an input Leaves the menu
	O/T / ins	Places insert marker (in text entry)
	CHANGE / clr	Deletes input (in text entry)
	PRINT	Inputs a point (in text entry)

Tab. 6.1  
Navigating through  
menu levels /  
Making an input

## 6.2

### Configuration menu

- Switch the balance off.
- First press the  key briefly and then immediately press the  key with a long press:

Menu items	Function
LANGUAGE	Selecting the language
UNIT-1	Unit in which the weighing results are displayed
SET DATA PRINT	Print formats, type of values to be printed out
SET CALIBRATION	Calibration method
SET WEIGHING MODE	Stability mode, Auto-Standby mode, zero correction, tare method, zero range
SET INTERFACE	Baud rate, parity, handshake functions of the peripheral interface
SET DATE AND TIME	Date and time
PASSWORD	Password protection for menu definitions
THEFTCODE	Activating / deactivating and changing the anti-theft code
KEY TONE	Keypad sound
BUS	Enable or disable the Precisa BUS or other accessories

Tab. 6.2  
Contents of the  
configuration menu

## 6.3

### Application menu

- Check that the balance is switched on.
- Press the  key with a long press:

Menu items	Function
SET APP. UNITS	Select an application
SETUP APPLICATION	Specify the parameters for the application selected under "APPLICATION"
SET STATISTICS	Statistics and storage functions
SET CHECK +/-	Define the nominal weight and limits for comparison weighings
AUTO-START ON	The selected application can, if required, be loaded automatically every time the balance is switched on

Tab. 6.3  
Contents of the  
application menu

## 7. Key menus

### 7.1 Tare menu

#### 7.1.1 Taring the balance

- Press the  key with a short press:  
→ The balance is tared.

#### 7.1.2 Other functions

- Press the  key until you reach the required function:

Tab. 7.1.2  
Contents of the  
tare menu

Menu items	Function
CALIBRATION	The balance begins the calibration process.

- Release the  key.

## 7.2

### Mode menu

#### 7.2.1

#### Switch to an activated application or back to weighing mode

- Press the  key with a short press:  
→ The balance switches to the last application activated. This acts like a switchover.

#### 7.2.2

#### Other functions

- Press the  key until you reach the required function:

Menu items	Function
WEIGHING	Switches to weighing mode
UNITS (COUNT; etc.)	Switches to the selected application units, counting, etc.
STATISTICS	Switches to the statistics program
CHECK +/-	Switches to the program for comparison weighings

- Release the  key.

Tab. 7.2.2  
Contents of the  
mode menu

### 7.3

## Print menu

#### 7.3.1

### Print individual value / report

- Press the  key with a short press:  
→ The individual value or report is printed out.

#### 7.3.2

### Other functions

- Press the  key until you reach the required function:

Menu items	Function
RESET PROD. COUNTER	Reset the product counter to 1
PRINT STATUS	Print out balance settings
PRINT APPLICATIONS	Print out application settings
PRINT CALIBRATIONS	Print out calibration information
PRINT FIRMWARE HIST	Print out firmware update information

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