



Moisture Applications
Ref No.: 090026

06.02.2009

Moisture determination of garden soil

This information is intended as a guideline for development of a moisture analysis method and parameter setting.

Materials vary throughout the world and we will be more than happy to test customers' individual samples.



Precisa Moisture Analysis XM 60 / XM 120 Halogen

Description of the sample

Utensils:
Preparation:
Sample weight (+/-10%):
Temperature program:
Temperature:
Stop-Criteria:
standby temperature

garden soil

aluminum pan, spoon
spread thinly on pan
2.0 g
standard
130 °C
Autostop 2 Digits [mg] / 30 s or 0.2 % / 30 s
80 °C

Moisture content:
Standard Deviation (+/-):
Average analysis time:

33.2% (average of 4 measurements)
0.8 %
7 min

Reference: