

## LAUNCH OF THE DRAFT SHIELD 321 LG

Starting from **January 1st 2018** the LG draft shield comes standard with all analytical balances of **series 320 XB<sup>BB</sup>, 321 LS and LT**. It is also compatible with series 321 LX, to be ordered as an option: 350-8723.


The new draft shield LG replaces the current draft shield XB at the same list price. Spare parts for the XB draft shield remain available as per usual.

### HIGHLIGHTS


- Improved metrological characteristics with better protection against air draft, especially when used in a fume hood
- Usability
  1. no stick-slip / jam effect
  2. better gliding characteristics
  3. improved quality
  4. improved visibility (no pillars in the front, glass in the back)
- Easy to clean
  1. side glasses are removeable without any tool and are dishwasher proof
  2. top glass is removeable without any tool and is dishwasher proof
  3. bottom plate is removeable and dishwasher proof
- Contemporary design



*Series 320 XB<sup>BB</sup>*



*Series 321 LS*



*Series 321 LT*