

OS UPGRADE FOR **prepASH®**

from Windows XP Embedded to Windows Embedded Standard 7

Windows 7 is the new standard operating system for the prepASH, reason being that Windows XP is no longer supported. As part of this upgrade the prepASH now comes with a new terminal as shown in the photo below.

The upgrade will automatically be applied to all purchase orders received from now on. Units purchased in the past can be upgraded to WES7 with the „**Upgrade kit 340 Terminal Windows 7**“, item no. **350-8899**.

All existing applications, settings, methods and results can easily be transferred to the new terminal.

WHAT IS NEW

- Operating System Windows from XPE to WES7
- New terminal
- Flat frame around the display, easier to clean
- Enhanced hardware
- New cover plate on the bottom side of the front panel

WHAT REMAINS THE SAME

- prepASH applications → fully compatible, no update necessary
- Connections and cables
- prepDATA and its functions
- Remote access and its functions
- Audit trail (since implementation of Software A1-1-2-05)

