


DRAFT SHIELD 321 LG-M, item no. 350-8731 suitable for 1mg, 0.01g, 0.1g and 1g balances

The LG draft shield launched for analytical balances of series 321 and 320 early this year is now available for other models (M/C/D/G) as well.

Suitable for

Series:	360 EP/ES, 321 LT/LX/LS and 320 XB ^{BB}
Capacity:	up to max. 2220g
Readability:	1mg, 0.01g, 0.1g, 1g



Ordering process

The purchase order should state the following item numbers depending on your choice of series and model:

- Series 360 and 321

Balances with SCS

1. xxx-xxxx-009 - e.g. EP 1220M: 360-9348-009 (instead of 360-9348-001)
2. 350-8731 Draft shield LG-M

Balances without SCS

1. xxx-xxxx-010 - e.g. LS 620C: 321-9459-010 (instead of 321-9459-002)
2. 350-8731 Draft shield LG-M

- LG draft shield for balances with standard item no. (incl. Series 320)

M-Balances

1. 350-8731 Draft shield LG-M
2. 350-7402 Adapter ring

C/D/G-Balances

1. 350-8731 Draft shield LG-M
2. 350-7402 Adapter ring
3. 320-7014 Weighing pan M 135 x 135 mm
4. 360-8014 Holder for weighing pan M

The applicable list prices for balances, the LG draft shield and all spare parts mentioned can be found in the current price lists. Your usual discounts will be applied.