


DRAFT SHIELD 321 LG-M suitable for 1mg, 0.01g, 0.1g and 1g balances

The LG draft shield launched for analytical balances of series 321 and 320 early this year is now available for other models (M/C/D/G) as well.

Suitable for

Series:	360 EP/ES, 321 LT/LX/LS and 320 XB ^{BB}
Capacity:	up to max. 2220g
Readability:	1mg, 0.01g, 0.1g, 1g



Ordering process

The purchase order should state the following item numbers depending on your choice of series and model:

- Series 360 and 321

Balances with SCS

1. xxx-xxxx-009 - e.g. EP 1220M: 360-9348-009 (instead of 360-9348-001)
2. 350-8731 Draft shield LG-M

Balances without SCS

1. xxx-xxxx-010 - e.g. LS 620C: 321-9459-010 (instead of 321-9459-002)
2. 350-8731 Draft shield LG-M

- Series 320
and

- Balances of Series 360/321 that were previously purchased and are already at your facility

M Balances

350-8732 Retrofit Kit Draft shield LG-M, M Balances

C/D/G Balances

350-8733 Retrofit Kit Draft shield LG-M, C/D/G Balances

The applicable list prices for balances, the LG draft shield and all spare parts mentioned can be found in the current price lists. Your usual discounts will be applied.