

New Printer LP 4024

Highlights

- Lightproof printout, accepted for use on official documents
- Single key operation
- Automatic interface settings detection
- Language diversity due to WPC1252 character set (Windows Latin-1)
- Two years warranty (Precisa standard)

Comprehensive scope of delivery

- Printer including international power supply
- Data cable DB9f / DB9m
- Data cable DB9f / RJ45
- 2 paper rolls
- 2 ribbons (1 assembled + 1 extra)
- Handbook brief



Specifications

Printing method	Impact Dot Matrix
Printer mechanism	EPSON M-192G
Printing characters per line	40* / 24 (5 x 7 dot matrix)
Character size	Approx. 1.7 mm (W) x 2.6 mm (H), 5 x 7 dot matrix
Printer paper	Normal paper, 57mm width, 50mm in diameter
Interface	D-SUB-9pin (RS-232C), male
Manipulation	1 key

*factory settings

Dimensions

