

Release Notes 1.2.0

New features

- Density application
- Add-on Checkweigher available for application Weighing, Counting, Percent and Dynamic weighing
- Wifi connection by USB Wifi adapter (350-8941)

Miscellaneous

- Soft button for levelling added
- Soft button for preset tare added
- Weight entries for the current application in the selected application unit
- Reference unit added to the System Settings menu for global weight entries
- Unit kilogram and pound added
- Info item Gross and Tare added
- Item for protocol format length 120 added
- Returning from Adjustment or Levelling to the main application by [Home] button
- Linearization with displayed default values

Bug fixes

- Unit Bath conversion factor corrected
- Unit Tola conversion factor corrected
- Unit Teal Singapore corrected
- Dynamic weighing tare/zero function corrected