

TEST CERTIFICATE

model: **LS 2200C**
S/N: 7204122

max. weight: 2200 g d: 0.01 g
e: 0.1 g

marked: yes no
int. calibration: yes no
accumulator: yes no

single range- dual range- floating range balance

corner load check	
nominal value = 800 g	limit: +/- 0.02 g
centre	0.00 g
front left	0.00 g
back left	0.00 g
back right	0.00 g
front right	0.00 g

calibration executed with weight set of	
calibration weight class	E2
serial number	76.03
approved until	2020-5

reproducing check	
series of 10 measurements	limit: 0.010 g
nominal value	2000 g
standard deviation	0.000 g

linearity check	
check value	limit: +/- 0.02 g
0 g	0.00 g
400 g	0.00 g
800 g	0.00 g
1200 g	0.00 g
1600 g	0.01 g
2000 g	0.00 g

measurements performed automatically

measurements performed automatically

Test passed

Date: 2020-02-17

Signature: MA  **Precisa Gravimetrics AG**
Inspected