

TEST CERTIFICATE

model: **HF 225SM-DR**
S/N: **8212345**

max. weight: 102 g / 225 g d: 0.00001 g / 0.0001 g
e: 0.001 g

marked: yes no
int. calibration: yes no
accumulator: yes no

single range- dual range- floating range balance

corner load check	
nominal value = 80 g	limit: 0.00020 g
centre	0.00000 g
front left	-0.00001 g
back left	0.00000 g
back right	0.00001 g
front right	0.00001 g

calibration executed with weight set of	
calibration weight class	E2
serial number	76.08
approved until	2019-6

reproducing check, fine	
series of 10 measurements	limit: 0.000030 g
nominal value	100 g
standard deviation	0.000011 g
measurements performed automatically	

linearity check, fine	
check value	limit: 0.00006 g
0 g	0.00000 g
40 g	0.00000 g
measurements performed automatically	

reproducing check, coarse	
series of 10 measurements	limit: 0.00010 g
nominal value	200 g
standard deviation	0.00001 g
measurements performed automatically	

linearity check, coarse	
check value	limit: 0.0002 g
0 g	0.0000 g
40 g	0.0000 g
80 g	0.0001 g
120 g	0.0000 g
160 g	-0.0001 g
200 g	0.0000 g
measurements performed automatically	

Test passed

Date: 2018-11-23

Signature: TV