



## INTRODUCING EEZYcount

Introducing the new EEZYcount tablet counter and weighing balance for all pharmacy, dispensary, veterinary and chemists premises in the UK and Ireland.

Conforming to the latest NAWI (EN45501/15) Directive for Non-Automated Weighing Instruments, the EEZYcount is trade stamped by our engineers in the United Kingdom. It is therefore suitable for sale and use in the United Kingdom & Republic of Ireland.

SWISS MADE  
SWISS TECHNOLOGY

## KEY FEATURES

- **Simple user operation:** seamlessly switch between counting and weighing mode with the press of one button
- **Internal calibration as standard:** provides assurance that instrument is in good working order, and ensures reliable and accurate results.
- **Cross contamination minimisation:** the use of removable, washable and disposable dispensing containers permits the operator to minimise the potential for cross contamination of different medicines
- **Balance protection:** with the plastic dust cover over the balance display & keyboard there is adequate protection from dust, spillages etc.
- **Stamped legal for trade:** to ensure compliance with EU Law

## BENEFITS

- **Increased accuracy due to the avoidance of human error:** the EEZYcount enables the correct quantity of pills or tablets to be measured.
- **Increased efficiency and cost savings:** the EEZYcount speeds up the process of counting, meaning that personnel time can be better utilised on other tasks.

**EEZYcount**  
SERIES 321 | LT 620C



# EEZYcount

SERIES 321 | LT 620C

## TECHNICAL SPECIFICATIONS

Model	Capacity	Readability	Repeatability	Linearity	Pan Size (mm)	CLASS / „e“	Part Number
LT 620C	620 g	0.01 g	0.005 g	10 mg	200 x 200	Class II / 0.1 g	340-9220-001

## SPARE PARTS

Part Number	Description
200-004	Scoops
350-8686	Dust Cover
320-7014	Pan
350-8384	Power Supply

## SERVICE OPTIONS

Description	Comments
Re-verify - Reverification of balance	
Service - Service and calibration of balance	Contact Precisa for price and information
Repair - In house repair of balance	

## CONTACT US

For price, service, calibration or any other queries please contact our sales office on +44 (0)1506 425300 from the UK and ROI

### Precisa Gravimetrics AG

Moosmattstrasse 32  
 8953 Dietikon, Switzerland  
 Tel. +41 44 744 28 28  
 Fax +41 44 744 28 38  
 Email: info@precisa.ch  
 www.precisa.com

