
Release Notes 1.0.6 (2026-02-18)

New features

- Introduced Session Log
- Added copy-and-paste functionality for result calculation formulas
- prepASH software version check (displays a hint when an update is needed)

Bug fixes

- Updated the result calculation mode identifier to replace rest with residue
- Fixed an issue where report creation was aborted if data was missing
- Resolved an issue with the calculation of method results
- Fixed result calculation when a comma is used as the decimal separator
- Fixed sorting of analyses by date
- Fixed display issues when screen scaling is set above 100%

Miscellaneous

- Compatible with PrepASH firmware version 1.3.0 or later

Release Notes 1.0.5

New features

- Support for prepSTATION 2.0
- Redesigned Sample List Editor

Bug fixes

- Corrected rounding errors in ASTM coal analysis results

Miscellaneous

- Compatible with PrepASH firmware version 1.3.0 or later

Release Notes 1.0.4

New features

- Support for ASTM coal analysis

Bug fixes

- Checkweigher information corrected
- Report: Only selected samples are displayed in the chart legend
- Report: Missing grid lines issue resolved
- Group name refresh issue fixed

Release Notes 1.0.3

Bug fixes

- Corrected calculation of additional results

Miscellaneous

- Compatible with PrepASH firmware version 1.0.4 or later

Release Notes 1.0.2

Miscellaneous

- Autostop (AS) information added to the measurement value table