



350-8719

Density Determination Kit, with Displacer,
390(SM/A)/520(A)

Dichtebestimmungssatz, mit Senkkörper,
390(SM/A)/520(A)



Scope of delivery – Lieferumfang



Pos.	Pcs.	Item number	Item	Artikel
1	1	350-7056	Thermometer	Thermometer
2	1	350-4024	Clip	Klammer
3	1	350-7055	Glass beaker	Becherglas
4	1	350-2066	Pedestal	Pult
5	1	350-7307	Support frame cpl.	Tragrahmen kpl.
6	1	350-7052	Double beaker	Doppelbecher
7	1	350-7194	Beaker with downholder, for samples with $\rho < 1\text{g/cm}^3$	Becher mit Niederhalter für Proben mit $\rho < 1\text{g/cm}^3$
8	1	350-7053	Beaker with hook	Becher mit Haken
9	1	350-7054	Displacer, $V=10\text{cm}^3$	Senkkörper, $V=10\text{cm}^3$